

小鼠抗 CBR1 单克隆抗体

中文名称： 小鼠抗 CBR1 单克隆抗体

英文名称： Anti-CBR1 mouse monoclonal antibody

别名： carbonyl reductase 1; CBR; hCBR1; SDR21C1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： CBR1

反应种属： Human

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	The protein encoded by this gene belongs to the short-chain dehydrogenases/reductases (SDR) family, which function as NADPH-dependent oxidoreductases having wide specificity for carbonyl compounds, such as quinones, prostaglandins, and various xenobiotics. Alternatively spliced transcript variants have been found for this gene.
Applications:	WB, IF
Name of antibody:	CBR1
Immunogen:	Fusion protein of human CBR1
Full name:	carbonyl reductase 1
Synonyms:	CBR; hCBR1; SDR21C1

SwissProt:	P16152
WB Predicted band size:	30 kDa
WB Positive control:	MDA-MB-468, A431 and Hela cell lysates
WB Recommended dilution:	1000-5000
IF positive control:	Hela cells
IF Recommend dilution:	50-100



