

小鼠抗 CBX8 单克隆抗体

中文名称: 小鼠抗 CBX8 单克隆抗体

英文名称: Anti-CBX8 mouse monoclonal antibody

别名: PC3; RC1

抗原: CBX8

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	Component of a Polycomb group (PcG) multiprotein PRC1-like complex, a complex class required to maintain the transcriptionally repressive state of many genes, including Hox genes, throughout development. PcG PRC1 complex acts via chromatin remodeling and modification of histones; it mediates monoubiquitination of histone H2A 'Lys-119', rendering chromatin heritably changed in its expressibility.
Applications:	WB, IHC
Name of antibody:	CBX8
Immunogen:	Fusion protein of human CBX8
Full name:	chromobox homolog 8 (CBX8)

Synonyms:	PC3; RC1
SwissProt:	Q9HC52
IHC positive control:	carcinoma of human thyroid tissue
IHC Recommend dilution:	25-100
WB Predicted band size:	43 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
WB Recommended dilution:	1000-5000