

兔抗 STIP1 多克隆抗体

中文名称：兔抗 STIP1 多克隆抗体

英文名称：Anti-STIP1 rabbit polyclonal antibody

别名：STIP1; HOP; IEF-SSP-3521; P60; STI1; STI1L

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：STIP1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	HOP, also known as stress-induced phospho protein 1 (STIP1), is a co-chaperone to the major heat shock proteins, Hsp70 and Hsp90, and appears in early receptor complexes. Through mutual binding to both Hsp70 and Hsp90, HOP functions as an adaptor that can integrate Hsp70 and Hsp90 interactions. HOP is an abundant and highly conserved protein which is composed of three tetratricopeptide repeat (TPR) domains (TPR1, TPR2a and TPR2b) and two DP repeat domains (DP1 and DP2), whose function has not been fully resolved.
Applications:	WB

Name of antibody:	STIP1
Immunogen:	Fusion protein of human STIP1
Full name:	stress-induced-phosphoprotein 1
Synonyms :	STIP1; HOP; IEF-SSP-3521; P60; STI1; STI1L
SwissProt:	P31948
WB Predicted band size:	63 kDa
WB Positive control:	Jurkat cells
WB Recommended dilution:	500-2000

