

兔抗 STK10 多克隆抗体

中文名称：兔抗 STK10 多克隆抗体

英文名称：Anti-STK10 rabbit polyclonal antibody

别名：LOK; PRO2729

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：STK10

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a member of the Ste20 family of serine/threonine protein kinases, and is similar to several known polo-like kinase kinases. The protein can associate with and phosphorylate polo-like kinase 1, and overexpression of a kinase-dead version of the protein interferes with normal cell cycle progression. The kinase can also negatively regulate interleukin 2 expression in T-cells via the mitogen activated protein kinase kinase 1 pathway.
Applications:	WB
Name of antibody:	STK10

Immunogen:	Synthesized peptide derived from internal of human S TK10.
Full name:	serine/threonine kinase 10
Synonyms :	LOK; PRO2729
SwissProt:	O94804
WB Predicted band size:	112 kDa
WB Positive control:	Hela cells and Huvec cells lysates
WB Recommended dilution:	500-3000

