

兔抗 STK16 多克隆抗体

中文名称：兔抗 STK16 多克隆抗体

英文名称： Anti-STK16 rabbit polyclonal antibody

别名： KRCT; MPSK; TSF1; PKL12

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： STK16

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	STK16 (serine/threonine kinase 16), also known as KRCT, MPSK, TSF1 or PKL12, is a 305 amino acid lipid-anchored membrane protein that belongs to the superfamily of serine/threonine protein kinase s. Expressed ubiquitously at low levels, STK16 is a protein kinase that can catalytically phosphorylate both serine and threonine resi
Applications:	ELISA, IHC
Name of antibody:	STK16
Immunogen:	Fusion protein of human STK16
Full name:	serine/threonine kinase 16

Synonyms :	KRCT; MPSK; TSF1; PKL12
SwissProt:	O75716
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	20-100

