

兔抗 STK24 多克隆抗体

中文名称：兔抗 STK24 多克隆抗体

英文名称：Anti-STK24 rabbit polyclonal antibody

别名：MST3; STK3; MST3B; STE20; HEL-S-95

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：STK24

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	This gene encodes a serine/threonine protein kinase that functions upstream of mitogen-activated protein kinase (MAPK) signaling. The encoded protein is cleaved into two chains by caspases; the N-terminal fragment (MST3/N) translocates to the nucleus and promotes programmed cell death. There is a pseudo gene for this gene on chromosome X. Alternative splicing results in multiple transcript variants.
Applications:	WB, IHC
Name of antibody:	STK24
Immunogen:	Synthesized peptide derived from internal of human

	STK24.
Full name:	serine/threonine kinase 24
Synonyms :	MST3; STK3; MST3B; STE20; HEL-S-95
SwissProt:	Q9Y6E0
IHC positive control:	Human colon carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	49 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000



