

## 小鼠抗 CDKN2B 单克隆抗体

中文名称: 小鼠抗 CDKN2B 单克隆抗体

英文名称: Anti-CDKN2B mouse monoclonal antibody

别名: CDK4I; INK4B; MTS2; P15; p15INK4b; TP15

储存: 冷冻 (-20℃) 避光

抗原: CDKN2B

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

**Background:**

This gene lies adjacent to the tumor suppressor gene CDKN2A in a region that is frequently mutated and deleted in a wide variety of tumors. This gene encodes a cyclin-dependent kinase inhibitor, which forms a complex with CDK4 or CDK6, and prevents the activation of the CDK kinases, thus the encoded protein functions as a cell growth regulator that controls cell cycle G1 progression. The expression of this gene was found to be dramatically induced by TGF beta, which suggested its role in the TGF beta induced growth inhibition. Two alternatively spliced transcript variants of this gene, which encode disti

	nct proteins, have been reported. [provided by RefSeq].
<b>Applications:</b>	WB, IHC, FC
<b>Name of antibody:</b>	CDKN2B
<b>Immunogen:</b>	Fusion protein of human CDKN2B
<b>Full name:</b>	cyclin-dependent kinase inhibitor 2B (p15, inhibits CDK4) (CDKN2B), transcript variant 1
<b>Synonyms:</b>	CDK4I; INK4B; MTS2; P15; p15INK4b; TP15
<b>SwissProt:</b>	P42772
<b>IHC positive control:</b>	human pancreas tissue and carcinoma of human thyroid tissue; adenocarcinoma of human endometrium tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	15 kDa
<b>WB Positive control:</b>	HepG2, Hela, Jurkat, MDCK, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000