

小鼠抗 CDX2 单克隆抗体

- 中文名称： 小鼠抗 CDX2 单克隆抗体
- 英文名称： Anti-CDX2 mouse monoclonal antibody
- 别名： caudal type homeobox 2; CDX3; CDX-3; CDX2/AS
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Mouse
- 抗原： CDX2
- 反应种属： Human, Monkey, Mouse
- 标记物： Unconjugate
- 克隆类型： Mouse Monoclonal

技术规格

Background:	This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis.
Applications:	WB, IP, IF
Name of antibody:	CDX2
Immunogen:	Fusion protein of human CDX2
Full name:	caudal type homeobox 2

Synonyms:	CDX3; CDX-3; CDX2/AS
SwissProt:	Q99626
WB Predicted band size:	34 kDa
WB Positive control:	HepG2, 3T3, Hela, SW480, C2C12, COS7 and Caco2 cell lysates
WB Recommended dilution:	1000-5000
IF positive control:	Hela cells
IF Recommend dilution:	50-100





