

小鼠抗 CEP63 单克隆抗体

中文名称： 小鼠抗 CEP63 单克隆抗体

英文名称： Anti-CEP63 mouse monoclonal antibody

别名： centrosomal protein 63; SCKL6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a protein with six coiled-coil domains. The protein is localized to the centrosome, a non-membraneous organelle that functions as the major microtubule-organizing center in animal cells. Several alternatively spliced transcript variants have been found, but their biological validity has not been determined.
Applications:	WB
Name of antibody:	CEP63
Immunogen:	Fusion protein of human CEP63
Full name:	centrosomal protein 63
Synonyms:	SCKL6
SwissProt:	Q96MT8
WB Predicted band size:	81 KD



WB Positive control:	293T, A549, MCF7, Hela cell lysates
WB Recommended dilution:	500-2000