

兔抗 STK32C 多克隆抗体

中文名称：兔抗 STK32C 多克隆抗体

英文名称：Anti-STK32C rabbit polyclonal antibody

别名：PKE; YANK3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：STK32C

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the serine/threonine protein kinase family. The specific function of this kinase is not known.
Applications:	ELISA, IHC
Name of antibody:	STK32C
Immunogen:	Synthetic peptide of human STK32C
Full name:	serine/threonine kinase 32C
Synonyms :	PKE; YANK3
SwissProt:	Q86UX6
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer

IHC Recommend dilution:

20-100



