

## 小鼠抗 CERS2 单克隆抗体

中文名称： 小鼠抗 CERS2 单克隆抗体

英文名称： Anti-CERS2 mouse monoclonal antibody

别名： ceramide synthase 2; L3; LASS2; SP260; TMSG1

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： CERS2

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein that has sequence similarity to yeast longevity assurance gene 1. Mutation or overexpression of the related gene in yeast has been shown to alter yeast lifespan. The human protein may play a role in the regulation of cell growth. Alternatively spliced transcript variants encoding the same protein have been described.
<b>Applications:</b>	WB, IHC
<b>Name of antibody:</b>	CERS2
<b>Immunogen:</b>	Fusion protein of human CERS2
<b>Full name:</b>	ceramide synthase 2
<b>Synonyms:</b>	L3; LASS2; SP260; TMSG1

<b>SwissProt:</b>	Q96G23
<b>IHC positive control:</b>	Human lung cancer and Human breast tissue
<b>IHC Recommend dilution:</b>	500-1000
<b>WB Predicted band size:</b>	45 KD
<b>WB Positive control:</b>	Jurkat, K562, HepG2, A549, Hela cell lysates
<b>WB Recommended dilution:</b>	500-2000