

小鼠抗 CLP1 单克隆抗体

- 中文名称： 小鼠抗 CLP1 单克隆抗体
- 英文名称： Anti-CLP1 mouse monoclonal antibody
- 别名： cleavage and polyadenylation factor I subunit 1; HEAB; hClp1
- 相关类别： 一抗
- 储存： 冷冻（-20℃）
- 宿主： Mouse
- 抗原： CLP1
- 反应种属： Human, Rat
- 标记物： Unconjugate
- 克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a member of the Clp1 family. The encoded protein is a multifunctional kinase which is a component of the tRNA splicing endonuclease complex and a component of the pre-mRNA cleavage complex II. This protein is implicated in tRNA, mRNA, and siRNA maturation. Mutations in this gene are associated with pontocerebellar hypoplasia type 10 (PCH10). Alternatively splice transcript variants encoding different isoforms have been found for this gene.
Applications:	WB
Name of antibody:	CLP1
Immunogen:	Fusion protein of human CLP1
Full name:	cleavage and polyadenylation factor I subunit 1
Synonyms:	HEAB; hClp1



SwissProt:	Q92989
WB Predicted band size:	48 KD
WB Positive control:	HepG2, Hela, A549, PC12, MCF7 cell lysates
WB Recommended dilution:	500-2000