

## 兔抗 STMN3 多克隆抗体

中文名称：兔抗 STMN3 多克隆抗体

英文名称： Anti-STMN3 rabbit polyclonal antibody

别名： SCLIP

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： STMN3

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a protein which is a member of the stathmin protein family. Members of this protein family form a complex with tubulins at a ratio of 2 tubulins for each stathmin protein. Microtubules require the ordered assembly of alpha- and beta-tubulins, and formation of a complex with stathmin disrupts microtubule formation and function. A pseudogene of this gene is located on chromosome 22. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	STMN3

<b>Immunogen:</b>	Synthesized peptide derived from internal of human STMN3.
<b>Full name:</b>	stathmin-like 3
<b>Synonyms :</b>	SCLIP
<b>SwissProt:</b>	Q9NZ72
<b>WB Predicted band size:</b>	21 kDa
<b>WB Positive control:</b>	293 cells lysate
<b>WB Recommended dilution:</b>	500-3000

