

小鼠抗 CLU 单克隆抗体

中文名称： 小鼠抗 CLU 单克隆抗体

英文名称： Anti-CLU mouse monoclonal antibody

别名： AAG4; APOJ; CLI; KUB1; MGC24903; SGP-2; SGP2; SP-40; TRPM-2; TRPM2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： CLU

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a secreted chaperone that can under some stress conditions also be found in the cell cytosol. It has been suggested to be involved in several basic biological events such as cell death, tumor progression, and neurodegenerative disorders. Alternate splicing results in both coding and non-coding variants.
Applications:	WB, IHC
Name of antibody:	CLU
Immunogen:	Fusion protein of human CLU
Full name:	clusterin (CLU), transcript variant 2

Synonyms:	AAG4; APOJ; CLI; KUB1; MGC24903; SGP-2; SGP2; SP-40 ; TRPM-2; TRPM2
SwissProt:	P10909
IHC positive control:	human pancreas tissue and human lymphoma tissue; h uman tonsil tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	52 kDa
WB Positive control:	SVT2, COS7 cell lysates
WB Recommended dilution:	200-1000