

## 小鼠抗 CMPK1 单克隆抗体

中文名称: 小鼠抗 CMPK1 单克隆抗体

英文名称: Anti-CMPK1 mouse monoclonal antibody

别名: CK; CMK; CMPK; UMK; UMP-CMPK; UMPK

储存: 冷冻 (-20℃) 避光

抗原: CMPK1

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes one of the enzymes required for cellular nucleic acid biosynthesis. This enzyme catalyzes the transfer of a phosphate group from ATP to CMP, UMP, or dCMP, to form the corresponding diphosphate nucleotide. Alternate splicing results in both coding and non-coding transcript variants. [provided by RefSeq, Feb 2012].
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	CMPK1
<b>Immunogen:</b>	Fusion protein of human CMPK1
<b>Full name:</b>	cytidine monophosphate (UMP-CMP) kinase 1, cytosolic (CMPK1), transcript variant 1
<b>Synonyms:</b>	CK; CMK; CMPK; UMK; UMP-CMPK; UMPK
<b>SwissProt:</b>	P30085

<b>IHC positive control:</b>	human pancreas tissue and human colon tissue; adenocarcinoma of human colon tissue and human kidney tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	22 kDa
<b>WB Positive control:</b>	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	1000-5000