

兔抗 STX2 多克隆抗体

中文名称：兔抗 STX2 多克隆抗体

英文名称：Anti-STX2 rabbit polyclonal antibody

别名：EPM; EPIM; STX2A; STX2B; STX2C

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：STX2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The product of this gene belongs to the syntaxin/pimorphin family of proteins. The syntaxins are a large protein family implicated in the targeting and fusion of intracellular transport vesicles. The product of this gene regulates epithelial-mesenchymal interactions and epithelial cell morphogenesis and activation. Alternatively spliced transcript variants encoding different isoforms have been identified.
Applications:	ELISA, WB, IHC
Name of antibody:	STX2
Immunogen:	Synthetic peptide of human STX2

Full name:	syntaxin 2
Synonyms :	EPM; EPIM; STX2A; STX2B; STX2C
SwissProt:	P32856
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human brain
IHC Recommend dilution:	25-100
WB Predicted band size:	33 kDa
WB Positive control:	A172 cells
WB Recommended dilution:	500-2000



