

兔抗 STX7 多克隆抗体

中文名称：兔抗 STX7 多克隆抗体

英文名称： Anti-STX7 rabbit ployclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： STX7

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	May be involved in protein trafficking from the plasma membrane to the early endosome (EE) as well as in homotypic fusion of endocytic organelles. Mediates the endocytic trafficking from early endosomes to late endosomes and lysosomes.
Applications:	ELISA, IHC
Name of antibody:	STX7
Immunogen:	Fusion protein of human STX7
Full name:	syntaxin 7
SwissProt:	O15400
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

