

## 兔抗 STX11 多克隆抗体

中文名称：兔抗 STX11 多克隆抗体

英文名称： Anti-STX11 rabbit ployclonal antibody

别名： FHL4; HLH4; HPLH4

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： STX11

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the syntaxin family. Syntaxins have been implicated in the targeting and fusion of intracellular transport vesicles. This family member may regulate protein transport among late endosomes and the trans-Golgi network. Mutations in this gene have been associated with familial hemophagocytic lymphohistiocytosis.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	STX11
<b>Immunogen:</b>	Fusion protein of human STX11
<b>Full name:</b>	syntaxin 11

<b>Synonyms :</b>	FHL4; HLH4; HPLH4
<b>SwissProt:</b>	O75558
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	33 kDa
<b>WB Positive control:</b>	Mouse lung tissue
<b>WB Recommended dilution:</b>	500-2000



