

## 小鼠抗 CTAG1B 单克隆抗体

- 中文名称: 小鼠抗 CTAG1B 单克隆抗体
- 英文名称: Anti-CTAG1B mouse monoclonal antibody
- 别名: CT6.1; CTAG; CTAG1; ESO1; LAGE-2; LAGE2B; NY-ESO-1
- 抗原: CTAG1B
- 储存: 冷冻 (-20℃) 避光
- 宿主: Mouse
- 反应种属: Human
- 相关类别: 一抗
- 标记物: Unconjugate
- 克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is an antigen that is overexpressed in many cancers but that is also expressed in normal testis. This gene is found in a duplicated region of the X-chromosome and therefore has a neighboring gene of identical sequence. [provided by RefSeq, Jan 2012].
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	CTAG1B
<b>Immunogen:</b>	Fusion protein of human CTAG1B
<b>Full name:</b>	cancer/testis antigen 1B (CTAG1B)
<b>Synonyms:</b>	CT6.1; CTAG; CTAG1; ESO1; LAGE-2; LAGE2B; NY-ESO-1

<b>SwissProt:</b>	P78358
<b>IHC positive control:</b>	adenocarcinoma of human endometrium tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	18 kDa
<b>WB Positive control:</b>	293T and MCF7 cell lysates
<b>WB Recommended dilution:</b>	200-1000