

兔抗 STXBP2 多克隆抗体

中文名称：兔抗 STXBP2 多克隆抗体

英文名称： Anti-STXBP2 rabbit polyclonal antibody

别名： FHL5; UNC18B; Hunc18b; UNC18-2; pp10122; MUNC18-2

抗原： STXBP2

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the STXBP/unc-18/SEC1 family. The encoded protein is involved in intracellular trafficking, control of SNARE (soluble NSF attachment protein receptor) complex assembly, and the release of cytotoxic granules by natural killer cells. Mutations in this gene are associated with familial hemophagocytic lymphohistiocytosis. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.
Applications:	ELISA, IHC

Name of antibody:	STXBP2
Immunogen:	Fusion protein of human STXBP2
Full name:	syntaxin binding protein 2
Synonyms :	FHL5; UNC18B; Hunc18b; UNC18-2; pp10122; MUNC18-2
SwissProt:	Q15833
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human liver cancer
IHC Recommend dilution:	25-100

