

兔抗 SULF1 多克隆抗体

中文名称：兔抗 SULF1 多克隆抗体

英文名称：Anti-SULF1 rabbit polyclonal antibody

别名：SULF-1; HSULF-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SULF1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an extracellular heparan sulfate endosulfatase. The encoded enzyme selectively removes 6-O-sulfate groups from heparan sulfate chains of heparan sulfate proteoglycans (HSPGs). The enzyme is secreted through the Golgi and is subsequently localized to the cell surface. The expression of this gene may be down-regulated in several types of cancer, including hepatocellular (HCC), ovarian and breast cancers. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC

Name of antibody:	SULF1
Immunogen:	Synthetic peptide of human SULF1
Full name:	sulfatase 1
Synonyms :	SULF-1; HSULF-1
SwissProt:	Q8IWU6
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	10-50

