

兔抗 SULF2 多克隆抗体

中文名称：兔抗 SULF2 多克隆抗体

英文名称： Anti-SULF2 rabbit polyclonal antibody

别名： HSULF-2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SULF2

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Heparan sulfate proteoglycans (HSPGs) act as coreceptors for numerous heparin-binding growth factors and cytokines and are involved in cell signaling. Heparan sulfate 6-O-endosulfatases, such as SULF2, selectively remove 6-O-sulfate groups from heparan sulfate. This activity modulates the effects of heparan sulfate by altering binding sites for signaling molecules (Dai et al., 2005 [PubMed 16192265]).
Applications:	ELISA, IHC
Name of antibody:	SULF2
Immunogen:	Synthetic peptide of human SULF2

Full name:	sulfatase 2
Synonyms:	HSULF-2
SwissProt:	Q8IWU5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	20-100

