

兔抗 SULT1B1 多克隆抗体

中文名称：兔抗 SULT1B1 多克隆抗体

英文名称： Anti-SULT1B1 rabbit polyclonal antibody

别名： ST1B1; ST1B2; SULT1B2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SULT1B1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Sulfotransferase enzymes catalyze the sulfate conjugation of many hormones, neurotransmitters, drugs, and xenobiotic compounds. These cytosolic enzymes are different in their tissue distributions and substrate specificities. The gene structure (number and length of exons) is similar among family members. However, the total genomic length of this gene is greater than that of other SULT1 genes.
Applications:	ELISA, WB, IHC
Name of antibody:	SULT1B1
Immunogen:	Fusion protein of human SULT1B1
Full name:	sulfotransferase family, cytosolic, 1B, member 1

Synonyms :	ST1B1; ST1B2; SULT1B2
SwissProt:	O43704
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human colon cancer and human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	35 kDa
WB Positive control:	Human fetal intestines tissue
WB Recommended dilution:	1000-5000



