

兔抗 SULT4A1 多克隆抗体

中文名称：兔抗 SULT4A1 多克隆抗体

英文名称：Anti-SULT4A1 rabbit polyclonal antibody

别名：sulfotransferase family 4A member 1; NST; BRSTL1; SULTX3; BR-STL-1; DJ388M5.3; hBR-STL-1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SULT4A1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the sulfotransferase family. The encoded protein is a brain-specific sulfotransferase believed to be involved in the metabolism of neurotransmitters. Polymorphisms in this gene may be associated with susceptibility to schizophrenia.
Applications:	ELISA, IHC
Name of antibody:	SULT4A1
Immunogen:	Fusion protein of human SULT4A1
Full name:	sulfotransferase family 4A member 1

Synonyms:	NST; BRSTL1; SULTX3; BR-STL-1; DJ388M5.3; hBR-STL-1
SwissProt:	Q9BR01
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	50-300

