

## 小鼠抗 DCP1B 单克隆抗体

中文名称： 小鼠抗 DCP1B 单克隆抗体

英文名称： Anti-DCP1B mouse monoclonal antibody

别名： decapping mRNA 1B; DCP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： DCP1B

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of a family of proteins that function in removing the 5' cap from mRNAs, which is a step in regulated mRNA decay. This protein localizes to cytoplasmic foci which are the site of mRNA breakdown and turnover. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	DCP1B
<b>Immunogen:</b>	Fusion protein of human DCP1B
<b>Full name:</b>	decapping mRNA 1B
<b>Synonyms:</b>	DCP1
<b>SwissProt:</b>	Q8IZD4
<b>WB Predicted band size:</b>	68 KD



<b>WB Positive control:</b>	Ramos, Hela, MCF7 cell lysates
<b>WB Recommended dilution:</b>	500-2000