

## 小鼠抗 DDT 单克隆抗体

中文名称: 小鼠抗 DDT 单克隆抗体

英文名称: Anti-DDT mouse monoclonal antibody

别名: DDCT

储存: 冷冻 (-20℃) 避光

抗原: DDT

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	D-dopachrome tautomerase converts D-dopachrome into 5,6-dihydroxyindole. The DDT gene is related to the migration inhibitory factor (MIF) in terms of sequence, enzyme activity, and gene structure. DDT and MIF are closely linked on chromosome 22. [provided by RefSeq, Jul 2008].
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	DDT
<b>Immunogen:</b>	Fusion protein of human DDT
<b>Full name:</b>	D-dopachrome tautomerase (DDT), transcript variant 1
<b>Synonyms:</b>	DDCT
<b>SwissProt:</b>	P30046

<b>IHC positive control:</b>	human prostate tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	13 kDa
<b>WB Positive control:</b>	HepG2, Hela, COS7, Jurkat, MDCK cell lysates
<b>WB Recommended dilution:</b>	500-2000