

## 兔抗 SUZ12 多克隆抗体

中文名称： 兔抗 SUZ12 多克隆抗体

英文名称： Anti-SUZ12 rabbit polyclonal antibody

别名： CHET9; JAZ1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SUZ12

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This zinc finger gene has been identified at the break points of a recurrent chromosomal translocation reported in endometrial stromal sarcoma. Recombination of these breakpoints results in the fusion of this gene and JAZF1. The protein encoded by this gene contains a zinc finger domain in the C terminus of the coding region.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SUZ12
<b>Immunogen:</b>	Synthetic peptide of human SUZ12
<b>Full name:</b>	SUZ12 polycomb repressive complex 2 subunit

<b>Synonyms:</b>	CHET9; JAZ1
<b>SwissProt:</b>	Q15022
<b>ELISA Recommended dilution:</b>	1000-5000
<b>WB Predicted band size:</b>	90 kDa
<b>WB Positive control:</b>	MCF-7, Hela and NCCIT cells
<b>WB Recommended dilution:</b>	1000-2000

