

小鼠抗 DEF6 单克隆抗体

中文名称: 小鼠抗 DEF6 单克隆抗体

英文名称: Anti-DEF6 mouse monoclonal antibody

别名: IBP; SLAT; SWAP70L

储存: 冷冻 (-20℃) 避光

抗原: DEF6

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	DEF6, or IBP, is a guanine nucleotide exchange factor (GEF) for RAC (MIM 602048) and CDC42 (MIM 116952) that is highly expressed in B and T cells (Gupta et al., 2003 [PubMed 12923183]).
Applications:	WB, IHC, IF
Name of antibody:	DEF6
Immunogen:	Fusion protein of human DEF6
Full name:	DEF6 guanine nucleotide exchange factor
Synonyms:	IBP; SLAT; SWAP70L
SwissProt:	Q9H4E7
IHC positive control:	human lymphoma tissue and human tonsil tissue and

	human bladder tissue
IHC Recommend dilution:	30-150
WB Predicted band size:	74 kDa
WB Positive control:	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
WB Recommended dilution:	500-2000
IF positive control:	COS7 cells transiently transfected by pCMV6-ENTRY
IF Recommend dilution:	100-200