

小鼠抗 DGKZ 单克隆抗体

中文名称： 小鼠抗 DGKZ 单克隆抗体

英文名称： Anti-DGKZ mouse monoclonal antibody

别名： diacylglycerol kinase zeta; DAGK5; DAGK6; DGK-ZETA; hDGKzeta

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： DGKZ

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene belongs to the eukaryotic diacylglycerol kinase family. It may attenuate protein kinase C activity by regulating diacylglycerol levels in intracellular signaling cascade and signal transduction. Alternative splicing occurs at this locus and multiple transcript variants encoding distinct isoforms have been identified.
Applications:	WB, IHC
Name of antibody:	DGKZ
Immunogen:	Fusion protein of human DGKZ
Full name:	diacylglycerol kinase zeta
Synonyms:	DAGK5; DAGK6; DGK-ZETA; hDGKzeta

SwissProt:	Q13574
IHC positive control:	Human lymphoma tissue
IHC Recommend dilution:	200-500
WB Predicted band size:	124 KD
WB Positive control:	Jurkat, Hela cell lysates
WB Recommended dilution :	200-1000