

兔抗 SYCP3 多克隆抗体

中文名称：兔抗 SYCP3 多克隆抗体

英文名称：Anti-SYCP3 rabbit polyclonal antibody

别名：COR1; SCP3; SPGF4

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SYCP3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes an essential structural component of the synaptonemal complex. This complex is involved in synapsis, recombination and segregation of meiotic chromosomes. Mutations in this gene are associated with azoospermia in males and susceptibility to pregnancy loss in females. Alternate splicing results in multiple transcript variants that encode the same protein.
Applications:	ELISA, WB, IHC
Name of antibody:	SYCP3
Immunogen:	Synthetic peptide of human SYCP3
Full name:	synaptonemal complex protein 3

Synonyms :	COR1; SCP3; SPGF4
SwissProt:	Q8IZU3
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	28 kDa
WB Positive control:	Jurkat cells and human seminoma tissue, human liver cancer tissue
WB Recommended dilution:	200-1000



