

## 小鼠抗 Digoxin 单克隆抗体

中文名称： 小鼠抗 Digoxin 单克隆抗体

英文名称： Anti-Digoxin mouse monoclonal antibody

别名： Lanoxin

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： Digoxin

反应种属： Digoxin

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	Digoxin is a beta blocker drug originally derived from the foxglove plant, Digitalis lanata. Digoxin is used primarily to improve the pumping ability of the heart in congestive heart failure (CHF), and treat problems such as high blood pressure. It is also used to help normalize some dysrhythmias (abnormal types of heartbeat).
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	Digoxin
<b>Immunogen:</b>	Digoxin-BSA
<b>Full name:</b>	Digoxin
<b>Synonyms:</b>	Lanoxin



<b>ELISA Recommended dilution:</b>	500-1000
<b>WB Recommended dilution:</b>	100-500