

## 小鼠抗 DIXDC1 单克隆抗体

中文名称: 小鼠抗 DIXDC1 单克隆抗体

英文名称: Anti-DIXDC1 mouse monoclonal antibody

别名: CCD1

抗原: DIXDC1

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human, Monkey

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a positive regulator of the Wnt signaling pathway. The encoded protein is found associated with gamma tubulin at the centrosome. Alternative splicing results in multiple transcript variants encoding different isoforms.
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	DIXDC1
<b>Immunogen:</b>	Fusion protein of human DIXDC1
<b>Full name:</b>	DIX domain containing 1 (DIXDC1), transcript variant 2
<b>Synonyms:</b>	CCD1
<b>SwissProt:</b>	Q155Q3

<b>IHC positive control:</b>	human breast tissue and adenocarcinoma of human breast tissue; human kidney tissue and human liver tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	77/53/24 kDa
<b>WB Positive control:</b>	HepG2, A549, COS7, Jurkat, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	200-1000