

小鼠抗 DNAJB6 单克隆抗体

中文名称： 小鼠抗 DNAJB6 单克隆抗体

英文名称： Anti-DNAJB6 mouse monoclonal antibody

别名： DnaJ heat shock protein family (Hsp40) member B6; DJ4; MRJ; DnaJ; HSJ2; HHDJ1; HSJ-2; MSJ-1; LGMD1D; LGMD1E

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Mouse

抗原： DNAJB6

反应种属： Human, Rat, Mouse

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	This gene encodes a member of the DNAJ protein family. DNAJ family members are characterized by a highly conserved amino acid stretch called the 'J-domain' and function as one of the two major classes of molecular chaperones involved in a wide range of cellular events, such as protein folding and oligomeric protein complex assembly. This family member may also play a role in polyglutamine aggregation in specific neurons. Alternative splicing of this gene results in multiple transcript variants; however, not all variants have been fully described.
Applications:	WB
Name of antibody:	DNAJB6

Immunogen:	Fusion protein of human DNAJB6
Full name:	DnaJ heat shock protein family (Hsp40) member B6
Synonyms:	DJ4; MRJ; DnaJ; HSJ2; HHDJ1; HSJ-2; MSJ-1; LGMD1D; LGM D1E
SwissProt:	O75190
WB Predicted band size:	36 KD
WB Positive control:	SH-SY5Y, HeLa, Jurkat, RAW264.7, PC12 cell lysates
WB Recommended dilution:	500-2000