

兔抗 SYNGR1 多克隆抗体

中文名称：兔抗 SYNGR1 多克隆抗体

英文名称： Anti-SYNGR1 rabbit polyclonal antibody

别名： synaptogyrin 1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SYNGR1

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes an integral membrane protein associated with presynaptic vesicles in neuronal cells. The exact function of this protein is unclear, but studies of a similar murine protein suggest that it functions in synaptic plasticity without being required for synaptic transmission. The gene product belongs to the synaptogyrin gene family. Three alternatively spliced variants encoding three different isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	SYNGR1

Immunogen:	Fusion protein of human SYNGR1
Full name:	synaptogyrin 1
SwissProt:	O43759
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	50-100

