

小鼠抗 DPPA4 单克隆抗体

中文名称: 小鼠抗 DPPA4 单克隆抗体

英文名称: Anti-DPPA4 mouse monoclonal antibody

别名: 2410091M23Rik

储存: 冷冻 (-20℃) 避光

抗原: DPPA4

宿主: Mouse

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	May be involved in the maintenance of active epigenetic status of target genes. May inhibit differentiation of embryonic cells into a primitive ectoderm lineage.
Applications:	WB, IF
Name of antibody:	DPPA4
Immunogen:	Fusion protein of human DPPA4
Full name:	developmental pluripotency associated 4 (DPPA4)
Synonyms:	2410091M23Rik
SwissProt:	Q7L190
WB Predicted band size:	34 kDa
WB Positive control:	HepG2, HeLa, SVT2, COS7, Jurkat, MDCK, MCF-7 cell lysate



	s
WB Recommended dilution:	500-2000