

兔抗 SYP 多克隆抗体

中文名称：兔抗 SYP 多克隆抗体

英文名称： Anti-SYP rabbit polyclonal antibody

别名： MRX96; MRXSYP

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： SYP

宿主： Rabbit

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes an integral membrane protein of small synaptic vesicles in brain and endocrine cells. The protein also binds cholesterol and is thought to direct targeting of vesicle-associated membrane protein 2 (synaptobrevin) to intracellular compartments. Mutations in this gene are associated with X-linked mental retardation (XLMR).
Applications:	WB, IHC
Name of antibody:	SYP
Immunogen:	Synthesized peptide derived from human Synaptophysin.
Full name:	synaptophysin

Synonyms :	MRX96; MRXSYP
SwissProt:	P08247
IHC positive control:	Human brain tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	34 kDa
WB Positive control:	LOVO cells lysate
WB Recommended dilution:	500-3000



