

## 小鼠抗 DSTN 单克隆抗体

中文名称： 小鼠抗 DSTN 单克隆抗体

英文名称： Anti-DSTN mouse monoclonal antibody

别名： ACTDP; ADF; bA462D18.2; HEL32

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： DSTN

反应种属： Human , Dog , Monkey

标记物： Unconjugate

克隆类型： mouse monoclonal

### 技术规格

<b>Background:</b>	The product of this gene belongs to the actin-binding proteins ADF family. This family of proteins is responsible for enhancing the turnover rate of actin in vivo. This gene encodes the actin depolymerizing protein that severs actin filaments (F-actin) and binds to actin monomers (G-actin). Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq].
<b>Applications:</b>	WB, IHC, IF, FC
<b>Name of antibody:</b>	DSTN
<b>Immunogen:</b>	Fusion protein of human DSTN

<b>Full name:</b>	dextrin (actin depolymerizing factor) (DSTN), transcript variant 1
<b>Synonyms:</b>	ACTDP; ADF; bA462D18.2; HEL32
<b>SwissProt:</b>	P60981
<b>IHC positive control:</b>	human pancreas tissue and carcinoma of human thyroid tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	19 kDa
<b>WB Positive control:</b>	HepG2, Hela, A549, COS7, Jurkat, MDCK, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000