

小鼠抗 DTNB 单克隆抗体

中文名称： 小鼠抗 DTNB 单克隆抗体

英文名称： Anti-DTNB mouse monoclonal antibody

别名： MGC17163; MGC57126

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Mouse

抗原： DTNB

反应种属： Human , Rat , Monkey , Mouse

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:

This gene encodes dystrobrevin beta, a component of the dystrophin-associated protein complex (DPC). The DPC consists of dystrophin and several integral and peripheral membrane proteins, including dystroglycans, sarcoglycans, syntrophins and dystrobrevin alpha and beta. The DPC localizes to the sarcolemma and its disruption is associated with various forms of muscular dystrophy. Dystrobrevin beta is thought to interact with syntrophin and the DP71 short form of dystrophin. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq].

Applications:	WB, IF, FC
Name of antibody:	DTNB
Immunogen:	Fusion protein of human DTNB
Full name:	dystrobrevin, beta (DTNB), transcript variant 5
Synonyms:	MGC17163; MGC57126
SwissProt:	O60941
WB Predicted band size:	71 kDa
WB Positive control:	HepG2, HeLa, SVT2, COS7, PC12, MCF-7 cell lysates
WB Recommended dilution:	200-1000