

## 兔抗 SYT6 多克隆抗体

中文名称：兔抗 SYT6 多克隆抗体

英文名称： Anti-SYT6 rabbit polyclonal antibody

别名： sytVI

抗原： SYT6

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene belongs to the synaptotagmin family. Synaptotagmins share a common domain structure that includes a transmembrane domain and a cytoplasmic region composed of 2 C2 domains, and are involved in calcium-dependent exocytosis of synaptic vesicles. This protein has been shown to be a key component of the secretory machinery involved in acrosomal exocytosis. Alternatively spliced transcript variants have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SYT6

<b>Immunogen:</b>	Fusion protein of human SYT6
<b>Full name:</b>	synaptotagmin VI
<b>Synonyms :</b>	sytVI
<b>SwissProt:</b>	Q5T7P8
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human ovarian cancer and human lung cancer
<b>IHC Recommend dilution:</b>	50-200

