

兔抗 SYTL4 多克隆抗体

中文名称：兔抗 SYTL4 多克隆抗体

英文名称： Anti-SYTL4 rabbit polyclonal antibody

别名： SLP4

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SYTL4

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the synaptotagmin like protein family. Members of this family are characterized by an N-terminal Rab27 binding domain and C-terminal tandem C2 domains. The encoded protein binds specific small Rab GTPases and is involved in intracellular membrane trafficking. This protein binds Rab27 and may be involved in inhibiting dense core vesicle exocytosis. Alternative splicing results in multiple transcript variants that encode the same protein.
Applications:	WB

Name of antibody:	SYTL4
Immunogen:	Synthesized peptide derived from internal of human SYTL4.
Full name:	synaptotagmin-like 4
Synonyms :	SLP4
SwissProt:	Q96C24
WB Predicted band size:	76 kDa
WB Positive control:	NIH/3T3 cells lysate
WB Recommended dilution:	500-3000

