

小鼠抗 EIF4E3 单克隆抗体

中文名称：小鼠抗 EIF4E3 单克隆抗体

英文名称： Anti-EIF4E3 mouse monoclonal antibody

别名： eukaryotic translation initiation factor 4E family member 3; eIF-4E3; eIF4E-3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： EIF4E3

反应种属： Human, Rat

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	EIF4E3 belongs to the EIF4E family of translational initiation factors that interact with the 5-prime cap structure of mRNA and recruit mRNA to the ribosome (Joshi et al., 2004 [PubMed 15153109]).[supplied by OMIM, Mar 20 08]
Applications:	WB
Name of antibody:	EIF4E3
Immunogen:	Fusion protein of human EIF4E3
Full name:	eukaryotic translation initiation factor 4E family member 3
Synonyms:	eIF-4E3; eIF4E-3
SwissProt:	Q8N5X7
WB Predicted band size:	24 KD



WB Positive control:	Hela, A549, PC12, MCF7 cell lysates
WB Recommended dilution:	200-1000