

兔抗 TAB3 多克隆抗体

中文名称：兔抗 TAB3 多克隆抗体

英文名称： Anti-TAB3 rabbit polyclonal antibody

别名： NAP1; MAP3K7IP3

相关类别： 一抗

抗原： TAB3

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

The product of this gene functions in the NF-kappa B signal transduction pathway. The encoded protein, and the similar and functionally redundant protein MAP3K7IP2/TAB2, forms a ternary complex with the protein kinase MAP3K7/TAK1 and either TRAF2 or T RAF6 in response to stimulation with the pro-inflammatory cytokines TNF or IL-1. Subsequent MAP3K7/TAK1 kinase activity triggers a signaling cascade leading to activation of the NF-kappaB transcription factor. The human genome contains a related pseudogene.

Applications:	ELISA, WB, IHC
Name of antibody:	TAB3
Immunogen:	Fusion protein of human TAB3
Full name:	TGF-beta activated kinase 1/MAP3K7 binding protein 3
Synonyms :	NAP1; MAP3K7IP3
SwissProt:	Q8N5C8
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human thyroid cancer and human brain
IHC Recommend dilution:	50-200
WB Predicted band size:	79 kDa
WB Positive control:	293T, Raw264.7, NIH/3T3, K562 and Jurkat cells
WB Recommended dilution:	1000-5000



