

小鼠抗 ELF1 单克隆抗体

中文名称：小鼠抗 ELF1 单克隆抗体

英文名称：Anti-ELF1 mouse monoclonal antibody

别名：E74 like ETS transcription factor 1; RIA1; EFTUD1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Mouse

抗原：ELF1

反应种属：Human

标记物：Unconjugate

克隆类型：mouse monoclonal

技术规格

Background:	This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	ELF1
Immunogen:	Fusion protein of human ELF1
Full name:	E74 like ETS transcription factor 1
Synonyms:	RIA1; EFTUD1



SwissProt:	P32519
WB Predicted band size:	67 KD
WB Positive control:	HUVEC cell lysate
WB Recommended dilution:	500-2000