

## 兔抗 TACO1 多克隆抗体

中文名称:	兔抗 TACO1 多克隆抗体
英文名称:	Anti-TACO1 rabbit polyclonal antibody
别名:	translational activator of cytochrome c oxidase I; CCDC44
储存:	冷冻 (-20℃)
抗原:	TACO1
宿主:	Rabbit
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
相关类别:	一抗
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a mitochondrial protein that functions as a translational activator of mitochondrially-encoded cytochrome c oxidase 1. Mutations in this gene are associated with Leigh syndrome.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TACO1
<b>Immunogen:</b>	Synthetic peptide of human TACO1
<b>Full name:</b>	translational activator of cytochrome c oxidase I
<b>Synonyms:</b>	CCDC44
<b>SwissProt:</b>	Q9BSH4
<b>ELISA Recommended dilution:</b>	500-1000

<b>IHC positive control:</b>	Human tonsil and Human thyroid cancer
<b>IHC Recommend dilution:</b>	20-100



