

小鼠抗 ENO2 单克隆抗体

中文名称： 小鼠抗 ENO2 单克隆抗体

英文名称： Anti-ENO2 mouse monoclonal antibody

别名： enolase 2; NSE; HEL-S-279

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： ENO2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： Mouse Monoclonal

技术规格

Background:	This gene encodes one of the three enolase isoenzymes found in mammals. This isoenzyme, a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development in rats and primates.
Applications:	WB
Name of antibody:	ENO2
Immunogen:	Fusion protein of human ENO2
Full name:	enolase 2
Synonyms:	NSE; HEL-S-279
SwissProt:	P09104

WB Predicted band size:	47 kDa
WB Positive control:	SHSY-5Y and U87-MG cell lysates
WB Recommended dilution:	1000-5000

