

兔抗 TAF4 多克隆抗体

中文名称：兔抗 TAF4 多克隆抗体

英文名称： Anti-TAF4 rabbit polyclonal antibody

别名： TAF2C; TAF4A; TAF2C1; TAFII130; TAFII135

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： TAF4

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Initiation of transcription by RNA polymerase II requires the activities of more than 70 polypeptides. The protein that coordinates these activities is transcription factor IID (TFIID), which binds to the core promoter to position the polymerase properly, serves as the scaffold for assembly of the remainder of the transcription complex, and acts as a channel for regulatory signals. TFIID is composed of the TATA-binding protein (TBP) and a group of evolutionarily conserved proteins known as TBP-associated factors or TAFs.
Applications:	WB

Name of antibody:	TAF4
Immunogen:	Synthesized peptide derived from internal of human TAF4.
Full name:	TAF4 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 135kDa
Synonyms :	TAF2C; TAF4A; TAF2C1; TAFII130; TAFII135
SwissProt:	O00268
WB Predicted band size:	110 kDa
WB Positive control:	NIH/3T3 cells lysate
WB Recommended dilution:	500-3000

