

兔抗 TAF5 多克隆抗体

中文名称：兔抗 TAF5 多克隆抗体

英文名称： Anti-TAF5 rabbit polyclonal antibody

别名： TAF2D; TAFII100

储存： 冷冻（-20℃） 避光

抗原： TAF5

宿主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Initiation of transcription by RNA polymerase II requires the activities of more than 70 polypeptides. The protein that coordinates these activities is transcription factor IID (TFIID), which binds to the core promoter to position the polymerase properly, serves as the scaffold for assembly of the remainder of the transcription complex, and acts as a channel for regulatory signals. TFIID is composed of the TATA-binding protein (TBP) and a group of evolutionarily conserved proteins known as TBP-associated factors or TAFs.
Applications:	WB, IHC

Name of antibody:	TAF5
Immunogen:	Synthesized peptide derived from internal of human TAF5.
Full name:	TAF5 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 100kDa
Synonyms :	TAF2D; TAFII100
SwissProt:	Q15542
IHC positive control:	Human liver carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	87 kDa
WB Positive control:	A549 cells lysate
WB Recommended dilution:	500-3000



