

## 兔抗 TAF12 多克隆抗体

- 中文名称:兔抗 TAF12 多克隆抗体
- 英文名称: Anti-TAF12 rabbit polyclonal antibody
- 别名: TAF2J; TAFII20
- 相关类别: 一抗
- 储存: 冷冻(-20℃)
- 宿 主: Rabbit
- 抗原: TAF12
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

## 技术规格

Background:	Control of transcription by RNA polymerase II invol ves the basal transcription machinery which is a col lection of proteins. These proteins with RNA polym erase II, assemble into complexes which are modul ated by transactivator proteins that bind to cis-reg ulatory elements located adjacent to the transcripti on start site. Some modulators interact directly wit h the basal complex, whereas others may act as bri dging proteins linking transactivators to the basal t
	ranscription factors. Some of these associated facto rs are weakly attached while others are tightly asso



	ciated with TBP in the TFIID complex. Among the l atter are the TAF proteins.
Applications:	ELISA, IHC
Name of antibody:	TAF12
Immunogen:	Fusion protein of human TAF12
Full name:	TAF12 RNA polymerase II, TATA box binding protei n (TBP)-associated factor, 20kDa
Synonyms :	TAF2J; TAFII20
SwissProt:	Q16514
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human esophagus cancer and human ovarian canc er
IHC Recommend dilution:	25-100

