

## 兔抗 TAF13 多克隆抗体

中文名称：兔抗 TAF13 多克隆抗体

英文名称： Anti-TAF13 rabbit polyclonal antibody

别名： TAF2K; TAFII18; TAFII-18; TAF(II)18

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： TAF13

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	Initiation of transcription by RNA polymerase II requires the activities of more than 70 polypeptides. The protein that coordinates these activities is transcription factor IID (TFIID), which binds to the core promoter to position the polymerase properly, serves as the scaffold for assembly of the remainder of the transcription complex, and acts as a channel for regulatory signals. TFIID is composed of the TATA-binding protein (TBP) and a group of evolutionarily conserved proteins known as TBP-associated factors or TAFs.
<b>Applications:</b>	IHC

<b>Name of antibody:</b>	TAF13
<b>Immunogen:</b>	Synthesized peptide derived from internal of human TAF13.
<b>Full name:</b>	TAF13 RNA polymerase II, TATA box binding protein (TBP)-associated factor, 18kDa
<b>Synonyms :</b>	TAF2K; TAFII18; TAFII-18; TAF(II)18
<b>SwissProt:</b>	Q15543
<b>IHC positive control:</b>	Human breast carcinoma tissue
<b>IHC Recommend dilution:</b>	50-100

