

兔抗 LRP1B 多克隆抗体

- 中文名称: 兔抗 LRP1B 多克隆抗体
- 英文名称: Anti-LRP1B rabbit polyclonal antibody
- 别名: LDL receptor related protein 1B; LRPDIT; LRP-DIT
- 相关类别: 一抗
- 储存: 冷冻 (-20℃)
- 宿主: Rabbit
- 抗原: LRP1B
- 反应种属: Human, Mouse
- 标记物: Unconjugate
- 克隆类型: rabbit polyclonal

技术规格

Background:	LRP1B belongs to the low density lipoprotein (LDL) receptor gene family. These receptors play a wide variety of roles in normal cell function and development due to their interactions with multiple ligands (Liu et al., 2001 [PubMed 11384978]).
Applications:	ELISA, IHC
Name of antibody:	LRP1B
Immunogen:	Synthetic peptide of human LRP1B
Full name:	LDL receptor related protein 1B
Synonyms:	LRPDIT; LRP-DIT

SwissProt:	Q9NZR2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

